

Производитель камеры	Модель камеры	SPECT или PLANAR	Тип источника	Размеры поверхности, в дюймах	Рекомендованная номинальная активность, мКи	Perflexion	Feather Lite	NES	C-Thru
Digirad	2020tc	SPECT	Square	8" x 8"	10	PF09R-057-10M	MED3722	N/A	CTRF10012
	Spect PAK	SPECT	Square	9" x 9"	10	PF09R-057-10M	MED3723	N/A	N/A
	Cardius / Ergo	SPECT	Square	8" x 6"	10	N/A	MED3724	N/A	CTRF10012
Dilon	6800	SPECT	Rectangular	8" x 6"	5	PF09R-057-5M	MED3725	N/A	CTRF10012
Elscont	SP4	SPECT	Circular	15,75"	10 max	PF18C-057-10M	MED3706	NES392	CTRF10050
	SP6	SPECT	Rectangular	21,25" x 15,75"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Helix	SPECT	Rectangular	21,25" x 15,75"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Cardial	SPECT	Rectangular	15,75" x 10"	10 max	PF16R-057-10M	MED3727	NES8510	CTRF10017
	Varicam	SPECT	Rectangular	21,25" x 15,75"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
GE Healthcare	300AM	PLANAR	Circular	10"	3	N/A	MED3700	NES298	N/A
	Maxi 2	PLANAR	Circular	15"	5	PF18C-057-5M	MED3706	NES391	CTRF10051
	Maxi 37	PLANAR	Circular	15"	5	PF18C-057-5M	MED3706	NES391	CTRF10051
	5000 Series	SPECT	Circular	24"	10	N/A	MED3708	NES8012	N/A
	Discovery VH	SPECT	Rectangular	20" x 14"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Discovery 530	SPECT	Rectangular	14,6" x 7,4"	10	N/A	MED3743	N/A	N/A
	Hawkeye	SPECT	Rectangular	20" x 14"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	DSTi/Dsi	SPECT	Rectangular	12,99" x 14"	10	N/A	MED3713	NES8480	N/A
	Maxxus	SPECT	Rectangular	21" x 16"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Starcam XRT	SPECT	Rectangular	21" x 16"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Starcam XCT	SPECT	Circular	15,4"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Starcam ACT	SPECT	Circular	15,4"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Optima	SPECT	Rectangular	14" x 9"	7,5	PF16R-057-7,5M	MED3727	NES8507	CTRF10017
	Millenium MPS	SPECT	Square	14" x 14"	10	N/A	MED3713	NES8480	N/A
	Millenium MPR, VG	SPECT	Rectangular	15,75" x 21,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Millenium MG	SPECT	Rectangular	14" x 20"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Neurocam	SPECT	Rectangular	7,9" x 6,7"	10	PF09R-057-10M	MED3722	N/A	CTRF10012
	Myosight	SPECT	Rectangular	20" x 14"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Infinia	SPECT	Rectangular	21,21" x 15,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Ventri	SPECT	Rectangular	14,6" x 7,4"	10	PF16R-057-10M	MED3727	NES8510	CTRF10017
IS²	SR (Single Head Rect.)	SPECT	Rectangular	20" x 14,8"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Pulse CDC/DCC	SPECT	Rectangular	15" x 10,2"	10	PF16R-057-10M	MED3727	NES8400	CTRF10017
	SC (Single Head Circ.)	SPECT	Circular	15"	10	PF18C-057-10M	MED3706	NES8510	CTRF10050
Mediso Medical	CardioSpect SC	SPECT	Rectangular	20,8" x 15,3"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	CardioSpect SR	SPECT	Rectangular	20,8" x 15,3"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	CardioSpect D90	SPECT	Rectangular	14,5" x 9"	10	N/A	N/A	N/A	CTRF13770 A
	CardioSpect VMAX	SPECT	Rectangular	20,8" x 15,3"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	NeuroCpect Quad	SPECT	Rectangular	9" x 8"	10	N/A	N/A	N/A	CTRF13770 A
	ARC 3000	SPECT	Circular	15"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Genesys	SPECT	Rectangular	20" x 15"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000

Philips Healthcare	Brightview	PLANAR	Rectangular	20,25" x 16"	10 max	PF24RW-057-10M	N/A	NES8610	CTRF10021
	Argus	SPECT	Rectangular	20" x 15"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Solus, Cardial, Vertex	SPECT	Rectangular	20" x 15"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Forte	SPECT	Rectangular	20" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Cirrus	SPECT	Circular	15"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Cardio MD	SPECT	Rectangular	9,2" x 15,4"	10 max	PF16R-057-10M	MED3727	NES8510	CTRF10017
	Cardio 60	SPECT	Rectangular	20" x 15"	10 max	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Precedence	SPECT/CT	Rectangular	20" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Skylight	SPECT	Rectangular	20" x 14"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
Philips Healthcare	Dynamo	PLANAR	Circular	10"	3	N/A	MED3700	NES297	N/A
	411	PLANAR	Circular	11"	5	N/A	MED3700	NES297	N/A
	412	PLANAR	Circular	12"	5	N/A	MED3700	NES297	N/A
	415	PLANAR	Circular	15"	5	PF18C-057-5M	MED3706	NES391	CTRF10051
	Prism 1000	SPECT	Rectangular	20" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Prism 2000	SPECT	Rectangular	20" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Prism 3000	SPECT	Rectangular	15,7" x 9,4"	10	PF16R-057-10M	MED3727	NES8510	CTRF10017
	SX300	SPECT	Rectangular	14" x 14"	10	N/A	MED3725	NES8480	N/A
	Axis/Irix	SPECT	Rectangular	21" x 15,5"	10	PF24R-057-10M	MED3709	NES8400	N/A
Siemens Medical Solutions	Multispect 2	SPECT	Square	21,25" x 15,5"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Multispect 3	SPECT	Rectangular	16" x 12"	10	N/A	MED3713	NES8480	N/A
	3700, 7500 Orbiter Series	SPECT	Circular	15,25"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Body Scan	SPECT	Rectangular	23,5" x 15,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Diacam	SPECT	Rectangular	21,25" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	3700, 7500 Orbiter Series	SPECT	Rectangular	21,25" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	C.cam	SPECT	Rectangular	14" x 8,4"	10	PF16R-057-10M	MED3727	NES8510	CTRF10017
	Symbia	SPECT	Rectangular	21,23" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	E.cam	SPECT	Rectangular	21,25" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
SMV (Sopha Medical)	DST	SPECT	Rectangular	11,7" x 15,8"	10	N/A	MED3713	NES8480	N/A
	DSX	SPECT	Rectangular	21,3" x 15,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	DS7	SPECT	Circular	15,75"	10	PF18C-057-10M	MED3706	NES392	CTRF10050
	Bodytrac	SPECT	Rectangular	21,3" x 15,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	DSTXL	SPECT	Rectangular	21,3" x 15,75"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	Vision FX Series	SPECT	Rectangular	20" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
Toshiba America	GCA7100A, 7200	SPECT	Rectangular	21,5" x 19"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	GCA 9300	SPECT	Rectangular	16,1" x 8,3"	10	PF16R-057-10M	MED3727	NES8510	CTRF10017
	GCA901, 901A	SPECT	Rectangular	20,1" x 14,5"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	GCA901WB	SPECT	Rectangular	20,1" x 14,6"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000
	GCA602, 602A	SPECT	Circular	13,75"	5	N/A	MED3700	NES297	N/A
	GCa601A	SPECT	Circular	13,75"	5	N/A	MED3700	NES297	N/A
	t.cam	SPECT	Rectangular	21,25" x 15"	10	PF24R-057-10M	MED3709	NES8400	CTRF10000